

GE VCT 64 SLICE PROTOCOLS

Protocol Name: CAROTID CTA	
INDICATIONS:	Evaluation of carotid arteries.
IV CONTRAST:	100 cc of nonionic contrast. 4cc/sec. 20ga or better IV.
ORAL CONTRAST:	None
PATIENT POSITIONING AND SCAN AREA:	Patient supine. Head in holder. Landmark at mid sternum. Cuts are taken from the mid aortic arch through the circle of willis. No gantry tilt.
GENERAL EXAM GUIDELINES:	2.5 x 2.5 mm axial images with 3 mm slab mip sag and cor recons. 250 DFOV.
TRIGGER OR CONTRAST DELAY:	Use smart prep or bolus timing or 15 second delay.
TECHNIQUE GUIDELINES:	120 kvp. Smart MA on.
LABEL:	SCOUT*, NECK*. The * means anything after these beginnings will hang correct in pacs
RECON 1:	.625mm thick/.625 mm interval. Std recon at 450ww/50wl. With 3 mm COR and SAG MPRs.
RECON 2:	2.5mm thick /2.5 mm interval. Std. recon 450ww/50wl.
ADDITIONAL POST PROCESSING:	3D carotids. Save 40 image rotation to pacs.
SEND EXAM:	Send scout, 2.5 axials, and SAG and COR slab MIP recons to pacs. Send 0.625 images to awwkstn1.
ARCHIVE:	
BILL AS:	Neck CTA

SCOUT PROTOCOL SET IN SCANNER

Scan Type	Start Loc	End Loc	kV	mA	Scout Plane	Voice Lights Timer
Scout	S250.00	I150.00	120	50 and 20	90 and 0	1
	Above circle of willis	Mid chest				

AXIAL IMAGE PROTOCOL SET IN SCANNER

Scan Type	Start Loc	End Loc	Thick Speed	Interval (mm)	Gantry Tilt	SFOV	kv	ma	DFOV	Show Recon1	Show Recon2	Show Recon3	Direct Vis
Helical Full 0.4 sec rotation	150.0	S200.0	0.625	0.625	S0.0	Small Body	120	550 12.60~	18.0	0.625mm Std	1.25m Std		On
	Mid chest	Above circle of willis	39.37 0.984:1					Smart mA On DR 40%		450ww 50wl For recons SS40	450ww 50w Send to Pacs SS40		3mm Cor and Sag Send to Pacs