

Protocol Name: **BREAST ROUTINE**

Body Part: **BREAST**

Coil: **INVIVO 16ch BREAST COIL**

Sequences	Seq Name	Image Plane	# Slices	TR	TE	Slice Thick	FOV	NSA	Matrix	Phase Dir
Localizer	gre	3 Plane	4	7.4	3.01	7	400	1	320	R>L
Axial T1	gre	Axial	176	6.1	2.6	1	340	1	384	R>L
Axial STIR	tse	Axial	55	5600	70	3	340	3	448	R>L
Cor STIR	tse	Coronal	48	4730	63	3	380	2	320	F>H
Sag T2 FS Right	tse	Sagittal	30	3010	64	4	180	2	320	H>F
Sag T2 FS Left	tse	Sagittal	30	3010	64	4	180	2	320	H>F
Dynamic 7Pass Post	fl3d_vibe	Axial	176	4.37	1.59	1	340	1	480	R>L

Details: Dynamic has 1 pre-gad pass, followed by 40 sec eval delay, then 7 passes following the injection