

Protocol Name: **PELVIS AND HIP** Body Part: **PELVIS / HIP** Coil: **BODY FLEX COIL**

Sequences	Seq Name	Image Plane	# Slices	TR	TE	Slice Thick	FOV	NSA	Matrix	Phase Dir
2 Plane Localizer	gre	2 Plane	10	7.8	3.69	6	400	1	256	R>L
Coronal T1	tse	Coronal	30	730	20	5	380	2	512	F>H
Coronal T2 FS	tse	Coronal	30	6470	53	5	380	1	512	F>H
Axial T1	tse	Axial	35	837	20	5	380	2	384	A>P
Axial T2 FS	tse	Axial	35	6230	51	5	380	1	384	A>P
Sag T2 FS	tse	Sagittal	25	3000	54	3.5	160	3	256	A>P
Small FOV side of interest										
Coronal T2 FS	tse	Coronal	24	5000	55	3	220	2	384	F>H
Axial T2 FS	tse	Axial	25	3000	54	3.5	180	3	256	A>P