

Protocol Name: **ELBOW ARTHROGRAM**

Body Part: **ELBOW**

Coil: **MEDIUM FLEX COIL**

Sequences	Seq Name	Image Plane	# Slices	TR	TE	Slice Thickness	FOV	NSA	Matrix	Phase Dir
Axial Localizer	gre	Axial	5	7.9	3.7	6	480	1	256	A>P
3 Plane Localizer	gre	3 Plane	9	8.3	3.7	6	200	1	256	F>H
Axial T2 FS	tse	Axial	30	3000	52	3	120	2	384	R>L
Axial T1	tse	Axial	30	500	12	3	120	2	384	R>L
Axial T1 FS	tse	Axial	30	500	11	3	120	2	320	R>L
Cor Obl T1 FS	tse	Coronal Obl	25	500	13	3	120	1	320	A>P
Cor Obl T2 FS	tse	Coronal Obl	25	3730	52	3	120	2	384	A>P
Sag Obl T1 FS	tse	Sagittal Obl	25	500	13	3	120	1	320	R>L
Sag Obl T2 FS	tse	Sagittal Obl	25	3730	52	3	120	2	384	R>L
Optional										
Cor T1 3D VIBE	fl3d_vibe	Coronal	144	11.3	4.9	0.5	120	1	288	A>P